


Product Datasheet

MRPL48 antibody (orb1313867)

Description	Mouse monoclonal antibody to MRPL48
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human MRPL48 (NP_057139) produced in E.coli.
Target	MRPL48
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	23.8 kDa
Uniprot ID	Q96GC5
NCBI	NM_016055
Dilution Range	WB 1:500~2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...