


Product Datasheet

BEGAIN antibody (orb1313878)

Description	Mouse monoclonal antibody to BEGAIN
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 374-517 of human BEGAIN (NP_001153003) produced in E.coli.
Target	BEGAIN
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9BUH8
NCBI	NM_001159531
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...