


## Product Datasheet

### Cyclin Y (CCNY) antibody (orb1313929)

|                            |                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to Cyclin Y (CCNY)                                      |
| <b>Species/Host</b>        | Mouse                                                                             |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                      |
| <b>Tested Applications</b> | WB                                                                                |
| <b>Immunogen</b>           | Full length human recombinant protein of human CCNY (NP_859049) produced in 293T. |
| <b>Target</b>              | CCNY                                                                              |
| <b>Note</b>                | For research use only                                                             |
| <b>Isotype</b>             | IgG2a                                                                             |
| <b>Clonality</b>           | Monoclonal                                                                        |
| <b>Uniprot ID</b>          | <a href="#">Q8ND76</a>                                                            |
| <b>NCBI</b>                | <a href="#">NM_181698</a>                                                         |
| <b>Dilution Range</b>      | WB 1:2000                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                   |



Western blot analysis of HEK293T cells w...