


Product Datasheet

THAP5 antibody (orb1313974)

Description	Mouse monoclonal antibody to THAP5
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human THAP5 (NP_872335) produced in E.coli.
Target	THAP5
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	26.4 kDa
Uniprot ID	Q7Z6K1
NCBI	NM_182529
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...