


Product Datasheet

TOMM40L antibody (orb1313983)

Description	Mouse monoclonal antibody to TOMM40L
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TOMM40L (NP_115550) produced in E.coli.
Target	TOMM40L
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	33.7 kDa
Uniprot ID	Q969M1
NCBI	NM_032174
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...