


Product Datasheet

DNAJB6 antibody (orb1313990)

Description	Mouse monoclonal antibody to DNAJB6
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human DNAJB6 (NP_005485) produced in HEK293T.
Target	DNAJB6
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	26.7 kDa
Uniprot ID	O75190
NCBI	NM_005494
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...