


Product Datasheet

VGLL1 antibody (orb1313991)

Description	Mouse monoclonal antibody to VGLL1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human VGLL1 (NP_057351) produced in E.coli.
Target	VGLL1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	28.5 kDa
Uniprot ID	Q99990
NCBI	NM_016267
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...