


Product Datasheet

Calcineurin A (PPP3CA) antibody (orb1314041)

Description	Mouse monoclonal antibody to Calcineurin A (PPP3CA)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 2-105 of human PPP3CA (NP_000935) produced in E.coli.
Target	PPP3CA
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q08209
NCBI	NM_000944
Dilution Range	WB 1:500, IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...