


Product Datasheet

CDX1 antibody (orb1314071)

Description	Mouse monoclonal antibody to CDX1
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CDX1 (NP_001795) produced in E.coli.
Target	CDX1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	28 kDa
Uniprot ID	P47902
NCBI	NM_001804
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...