


Product Datasheet

RPAIN antibody (orb1314112)

Description	Mouse monoclonal antibody to RPAIN
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human RPAIN (NP_001028174) produced in HEK293T cell.
Target	RPAIN
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q86UA6
NCBI	NM_001160243
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...