


Product Datasheet

MAP126 (SPAG5) antibody (orb1314151)

Description	Mouse monoclonal antibody to MAP126 (SPAG5)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-297 of human SPAG5 (NP_006452) produced in E.coli.
Target	SPAG5
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	134.2 kDa
Uniprot ID	Q96R06
NCBI	NM_006461
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...