


Product Datasheet

U2AF1L3 (U2AF1L4) antibody (orb1314160)

Description	Mouse monoclonal antibody to U2AF1L3 (U2AF1L4)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human U2AF1L4 (NP_659424) produced in E.coli.
Target	U2AF1L4
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	21.8 kDa
Uniprot ID	Q8WU68
NCBI	NM_144987
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...