


Product Datasheet

RNF149 antibody (orb1314168)

Description	Mouse monoclonal antibody to RNF149
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 250-383 of human RNF149 (NP_775918) produced in E.coli.
Target	RNF149
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	43 kDa
Uniprot ID	Q8NC42
NCBI	NM_173647
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...