


Product Datasheet

PGBD3 antibody (orb1314231)

Description	Mouse monoclonal antibody to PGBD3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 313-593 of human PGBD3 (NP_736609) produced in E.coli.
Target	PGBD3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q8N328
NCBI	NM_170753
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...