


Product Datasheet

SLC2A4RG antibody (orb1314240)

Description	Mouse monoclonal antibody to SLC2A4RG
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-259 of human SLC2A4RG (NP_064446) produced in E.coli.
Target	SLC2A4RG
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	41.1 kDa
Uniprot ID	Q9NR83
NCBI	NM_020062
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...