


Product Datasheet

LZTS1 antibody (orb1314272)

Description	Mouse monoclonal antibody to LZTS1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human LZTS1 (NP_066300) produced in HEK293T Cell.
Target	LZTS1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	66.4 kDa
Uniprot ID	Q9Y250
NCBI	NM_021020
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...