


Product Datasheet

ZFP37 antibody (orb1314282)

Description	Mouse monoclonal antibody to ZFP37
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ZFP37(NP_003399) produced in E.coli.
Target	ZFP37
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	71 kDa
Uniprot ID	Q9Y6Q3
NCBI	NM_003408
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...