


Product Datasheet

TSEN34 antibody (orb1314329)

Description	Mouse monoclonal antibody to TSEN34
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TSEN34(NP_076980) produced in E.coli.
Target	TSEN34
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	33.5 kDa
Uniprot ID	Q9BSV6
NCBI	NM_024075
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...