


Product Datasheet

ZNF18 antibody (orb1314363)

Description	Mouse monoclonal antibody to ZNF18
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 183-405 of human ZNF18(NP_653281) produced in E.coli.
Target	ZNF18
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	62.1 kDa
Uniprot ID	P17022
NCBI	NM_144680
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...