


Product Datasheet

CD299 (CLEC4M) antibody (orb1314396)

Description	Mouse monoclonal antibody to CD299 (CLEC4M)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 74-329 of human CLEC4M (NP_055072) produced in E.coli.
Target	CLEC4M
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	45.2 kDa
Uniprot ID	Q9H2X3
NCBI	NM_014257
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...