


Product Datasheet

ZNF543 antibody (orb1314405)

Description	Mouse monoclonal antibody to ZNF543
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ZNF543 (NP_998763) produced in E.coli.
Target	ZNF543
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	68.2 kDa
Uniprot ID	Q08ER8
NCBI	NM_213598
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...