

Product Datasheet

ST2 (IL1RL1) antibody (orb1314414)

Description	Mouse monoclonal antibody to ST2 (IL1RL1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Full-length human recombinant protein of human ST2(NP_057316) produced in HEK293E
Target	IL1RL1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	63.36 kDa
Uniprot ID	Q01638
NCBI	NM_016232
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.