


Product Datasheet

ESE1 (ELF3) antibody (orb1314430)

Description	Mouse monoclonal antibody to ESE1 (ELF3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 175-286 of human ELF3(NP_001107781) produced in E.coli.
Target	ELF3
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	41.3 kDa
Uniprot ID	P78545
NCBI	NM_001114309
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...