


Product Datasheet

TEKT2 antibody (orb1314452)

Description	Mouse monoclonal antibody to TEKT2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-244 of human TEKT2(NP_055281) produced in E.coli.
Target	TEKT2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	49.5 kDa
Uniprot ID	Q9UIF3
NCBI	NM_014466
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...