


Product Datasheet

AMDHD1 antibody (orb1314462)

Description	Mouse monoclonal antibody to AMDHD1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human AMDHD1 (NP_689648) produced in HEK293T cell.
Target	AMDHD1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	46.6 kDa
Uniprot ID	Q96NU7
NCBI	NM_152435
Dilution Range	WB 1:2000, IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...