


Product Datasheet

IQGAP1 antibody (orb1314483)

Description	Mouse monoclonal antibody to IQGAP1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 464-609 of human IQGAP1(NP_003861) produced in E.coli.
Target	IQGAP1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P46940
NCBI	NM_003870
Dilution Range	WB 1:500, IHC 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...