


Product Datasheet

MAG1 (GPAT3) antibody (orb1314565)

Description	Mouse monoclonal antibody to MAG1 (GPAT3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 41-131 of human AGPAT9 (NP_116106) produced in E.coli.
Target	GPAT3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q53EU6
NCBI	NM_032717
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...