


Product Datasheet

NUDT4 antibody (orb1314625)

Description	Mouse monoclonal antibody to NUDT4
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human NUDT4 (NP_061967) produced in E.coli.
Target	NUDT4
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	20.1 kDa
Uniprot ID	A0A024RBG1
NCBI	NM_019094
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...