

Product Datasheet

TAK1 (MAP3K7) antibody (orb1314654)

Description	Mouse monoclonal antibody to TAK1 (MAP3K7)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 254-367 of human MAP3K7(NP_003179) produced in E.coli.
Target	MAP3K7
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	64 kDa
Uniprot ID	O43318
NCBI	NM_003188
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

