


Product Datasheet

TRMT12 antibody (orb1314711)

Description	Mouse monoclonal antibody to TRMT12
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TRMT12 (NP_060426) produced in E.coli.
Target	TRMT12
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	50.1 kDa
Uniprot ID	Q53H54
NCBI	NM_017956
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...