


Product Datasheet

ZRANB2 antibody (orb1314723)

Description	Mouse monoclonal antibody to ZRANB2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 95-320 of human ZRANB2 (NP_005446) produced in E.coli.
Target	ZRANB2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	36.1 kDa
Uniprot ID	O95218
NCBI	NM_005455
Dilution Range	IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...