


Product Datasheet

ZCCHC17 antibody (orb1314732)

Description	Mouse monoclonal antibody to ZCCHC17
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ZCCHC17 (NP_057589) produced in E.coli.
Target	ZCCHC17
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	27.4 kDa
Uniprot ID	Q9NP64
NCBI	NM_016505
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...