


Product Datasheet

HOXD4 antibody (orb1314750)

Description	Mouse monoclonal antibody to HOXD4
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human HOXD4 (NP_055436) produced in E.coli.
Target	HOXD4
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	27.7 kDa
Uniprot ID	P09016
NCBI	NM_014621
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...