


Product Datasheet

RAB3A antibody (orb1314752)

Description	Mouse monoclonal antibody to RAB3A
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human RAB3A(NP_002857) produced in E.coli.
Target	RAB3A
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P20336
NCBI	NM_002866
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...