


Product Datasheet

C14orf50 (PPP1R36) antibody (orb1314769)

Description	Mouse monoclonal antibody to C14orf50 (PPP1R36)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human PPP1R36(NP_758953) produced in E.coli.
Target	PPP1R36
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	49.3 kDa
Uniprot ID	Q96LQ0
NCBI	NM_172365
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...