

Product Datasheet

MRPL42 antibody (orb1314822)

Description	Mouse monoclonal antibody to MRPL42
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human MRPL42(NP_054769) produced in E.coli.
Target	MRPL42
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	16.5 kDa
Uniprot ID	Q9Y6G3
NCBI	NM_014050
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...