


Product Datasheet

PITPNB antibody (orb1314842)

Description	Mouse monoclonal antibody to PITPNB
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human PITPNB(NP_036531) produced in E.coli.
Target	PITPNB
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	31.4 kDa
Uniprot ID	P48739
NCBI	NM_012399
Dilution Range	WB 1:500~2000, IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...