


Product Datasheet

C1orf83 (TCEANC2) antibody (orb1314854)

Description	Mouse monoclonal antibody to C1orf83 (TCEANC2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TCEANC2(NP_694580) produced in E.coli.
Target	TCEANC2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	24 kDa
Uniprot ID	Q96MN5
NCBI	NM_153035
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...