


Product Datasheet

TBX20 antibody (orb1314857)

Description	Mouse monoclonal antibody to TBX20
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-340 of human TBX20(NP_001071121) produced in E.coli.
Target	TBX20
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	49.1 kDa
Uniprot ID	Q9UMR3
NCBI	NM_001077653
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...