


Product Datasheet

ARHGAP12 antibody (orb1314904)

Description	Mouse monoclonal antibody to ARHGAP12
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 286-452 of human ARHGAP12 (NP_060757) produced in E.coli.
Target	ARHGAP12
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	96.1 kDa
Uniprot ID	Q8IWW6
NCBI	NM_018287
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...