


Product Datasheet

RIAM (APBB1IP) antibody (orb1314918)

Description	Mouse monoclonal antibody to RIAM (APBB1IP)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human APBB1IP (NP_061916) produced in HEK293T cell.
Target	APBB1IP
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	73 kDa
Uniprot ID	Q7Z5R6
NCBI	NM_019043
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...