


Product Datasheet

Nestin (NES) antibody (orb1314937)

Description	Mouse monoclonal antibody to Nestin (NES)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 99-349 of human NES(NP_006608) produced in E.coli.
Target	NES
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P48681
NCBI	NM_006617
Dilution Range	WB 1:2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...