

Product Datasheet

Estrogen Related Receptor alpha (ESRRA) antibody (orb1314943)

Description Mouse monoclonal antibody to Estrogen Related Receptor alpha (ESRRA)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Human recombinant protein fragment corresponding to amino acids 206-417 of human ESRRA (NP_004442) produced in E.coli.

Target ESRRA

Note For research use only

Isotype IgG2a

Clonality Monoclonal

MW 46 kDa

Uniprot ID [P11474](#)

NCBI [NM_004451](#)

Dilution Range WB 1:500~2000

Expiration Date 12 months from date of receipt.

