


Product Datasheet

DACT1 antibody (orb1314959)

Description	Mouse monoclonal antibody to DACT1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 289-555 of human DACT1(NP_057735) produced in E.coli.
Target	DACT1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9NYF0
NCBI	NM_016651
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...