


Product Datasheet

Vitamin D Binding protein (GC) antibody (orb1314963)

Description	Mouse monoclonal antibody to Vitamin D Binding protein (GC)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human GC (NP_000574) produced in HEK293T cell.
Target	GC
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	51.1 kDa
Uniprot ID	P02774
NCBI	NM_000583
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...