


Product Datasheet

ELP4 antibody (orb1314977)

Description	Mouse monoclonal antibody to ELP4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 133-407 of human ELP4(NP_061913) produced in E.coli.
Target	ELP4
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	47 kDa
Uniprot ID	Q96EB1
NCBI	NM_019040
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...