


Product Datasheet

MGMT antibody (orb1314993)

Description	Mouse monoclonal antibody to MGMT
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human MGMT (NP_002403) produced in HEK293T cell.
Target	MGMT
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P16455
NCBI	NM_002412
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...