


Product Datasheet

ZNF230 antibody (orb1315015)

Description	Mouse monoclonal antibody to ZNF230
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ZNF230 (NP_006291) produced in E.coli.
Target	ZNF230
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	54.4 kDa
Uniprot ID	Q9UIE0
NCBI	NM_006300
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...