


Product Datasheet

ZNF217 antibody (orb1315043)

Description	Mouse monoclonal antibody to ZNF217
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 730-1048 of human ZNF217(NP_006517) produced in E.coli.
Target	ZNF217
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	115.1 kDa
Uniprot ID	O75362
NCBI	NM_006526
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...