


Product Datasheet

COPS3 antibody (orb1315086)

Description	Mouse monoclonal antibody to COPS3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 224-423 of human COPS3(NP_003644) produced in E.coli.
Target	COPS3
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q9UNS2
NCBI	NM_003653
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...